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- 1). FDA-NIH BIOMARKER WORKING GROUP. BEST (Biomarkers, Endpoints, and other Tools) resource [Internet]. 2016.
- 2). Lee JW, Devanarayan V, Barrett YC, Weiner R, Allinson J, Fountain S, Keller S, Weinryb I, Green M, Duan L, Rogers JA, Millham R, O'Brien PJ, Sailstad J, Khan M, Ray C, Wagner JA. Fit-for-purpose method development and validation for successful biomarker measurement. *Pharm Res*. 2006;23(2):312-28.
- 3). Piccoli SP, Garofolo F. Biomarker assay validation. *Bioanalysis*. 2018;10(12):889-891.
- 4). Biomarker Assay Collaborative Evidentiary Considerations Writing Group, Critical Path Institute. Points to Consider Document: Scientific and Regulatory Considerations for the Analytical Validation of Assays Used in the Qualification of Biomarkers in Biological Matrices. June 2019. <https://c-path.org/wp-content/uploads/2019/06/evidconsid-whitepaper-analyticalsectionv2019.pdf>
- 5). Ohtsu Y, Tanaka S, Igarashi H, Kakehi M, Mori T, Nakamura T, Ohashi R, Shimizu H, Yasuda Y, Okayama T, Kakuo H, Yokoi H, Horiuchi M, Katashima M, Nakamura R, Saito K, Saito Y. Analytical method validation for biomarkers as a drug development tool: points to consider. *Bioanalysis*. 2021 Sep;13(18):1379-1389. doi: 10.4155/bio-2021-0173.
- 6). US FDA. Table of Surrogate Endpoints That Were the Basis of Drug Approval or Licensure. <https://www.fda.gov/drugs/development-resources/table-surrogate-endpoints-were-basis-drug-approval-or-licensure>
- 7). Schofield AL, Brown JP, Brown J, Wilczynska A, Bell C, Glaab WE, Hackl M, Howell L, Lee S, Dear JW, Remes M, Reeves P, Zhang E, Allmer J, Norris A, Falciani F, Takeshita LY, Seyed Forootan S, Sutton R, Park BK, Goldring C. Systems analysis of miRNA biomarkers to inform drug safety. *Arch Toxicol*. 2021;95(11):3475-3495.
- 8). Goossens N, Nakagawa S, Sun X, Hoshida Y. Cancer biomarker discovery and validation. *Transl Cancer Res*. 2015;4(3):256-269.
- 9). ICH. Nonclinical Biodistribution Considerations For Gene Therapy Products S12. 2023. [https://database.ich.org/sites/default/files/ICH\\_S12\\_Step4\\_Guideline\\_2023\\_0314.pdf](https://database.ich.org/sites/default/files/ICH_S12_Step4_Guideline_2023_0314.pdf)
- 10). ICH. General Principles to Address Virus and Vector Shedding. [https://admin.ich.org/sites/default/files/2019-04/ICH\\_Considerations\\_Viral-Vector\\_Shedding\\_.pdf](https://admin.ich.org/sites/default/files/2019-04/ICH_Considerations_Viral-Vector_Shedding_.pdf)
- 11). Bustin SA, Benes V, Garson JA, Hellemans J, Huggett J, Kubista M, Mueller R, Nolan T, Pfaffl MW, Shipley GL, Vandesompele J, Wittwer CT. The MIQE guidelines: minimum information for publication of quantitative real-time PCR experiments. *Clin Chem*. 2009;55(4):611-22.
- 12). Hays A, Islam R, Matys K, Williams D. Best Practices in qPCR and dPCR Validation in Regulated Bioanalytical Laboratories. *AAPS J*. 2022;24(2):36.
- 13). Corsaro B, Yang TY, Murphy R, Sonderegger I, Exley A, Bertholet S, Dakappagari N, Dessy F, Garofolo F, Kierstead L, Koch H, Sarikonda G, Savoie N, Siggers R, Solstad T, Lu Y, Milton M, Marshall JC, DelCarpini J, Gorovits B, Gupta S, Jesaitis L, Kamerud J, Kromminga A, Ma A, McNally J, Yan H, Wu B, Verthelyi D, Kirshner S, Pedras-Vasconcelos J, Rajadhyaksha M, Staack RF, Cherry E, Cludts I, Dahlbäck M, Gunn GR, Ishii-Watabe A, Jawa V, Kubiak R, Partridge M,

- Petrillo M, Pine SO, Poetzl J, Song S, Stebbins C, Wu Y, Zhang L, Kar S, Liang M, Abhari MR, Schweighardt B, Stubenrauch K, Xu Y. 2020 White Paper on Recent Issues in Bioanalysis: Vaccine Assay Validation, qPCR Assay Validation, QC for CAR-T Flow Cytometry, NAb Assay Harmonization and ELISpot Validation (Part 3 - Recommendations on Immunogenicity Assay Strategies, NAb Assays, Biosimilars and FDA/EMA Immunogenicity Guidance/Guideline, Gene & Cell Therapy and Vaccine Assays). *Bioanalysis*. 2021;13(6):415-463.
- 14). Wissel M, Poirier M, Satterwhite C, Lin J, Islam R, Zimmer J, Khadang A, Zemo J, Lester T, Fjording M, Hays A, Hughes N, Garofolo F, Guilbaud R, Groeber E, Renfrew H, Colletti K, Yu M, Lin J, Fang X, Shah S, Garofolo W, Kar S, Hayes R, Pirro J, Kane C, Luna M, Xu A, Cape S, O'Dell M, Wheller R, Ritzen H, Vance J, Farley E, Matys K, Tabler E, Mylott W, Yuan M, Karnik S, Voelker T, DuBey I, Williard C, Shi J, Yamashita J. Recommendations on qPCR/ddPCR assay validation by GCC. *Bioanalysis*. 2022;14(12):853-863.
  - 15). Ma H, Bell KN, Loker RN. qPCR and qRT-PCR analysis: Regulatory points to consider when conducting biodistribution and vector shedding studies. *Mol Ther Methods Clin Dev*. 2020;20:152-168.
  - 16). Uchiyama A, Naritomi Y, Hashimoto Y, Hanada T, Watanabe K, Kitta K, Suzuki G, Komatsuno T, Nakamura T. Understanding quantitative polymerase chain reaction bioanalysis issues before validation planning: Japan Bioanalysis Forum discussion group. *Bioanalysis*, 2023 Jan 9. doi: 10.4155/bio-2022-0190.
  - 17). US FDA. Guidelines for the Validation of Analytical Methods for Nucleic Acid Sequence-Based Analysis of Food, Feed, Cosmetics and Veterinary Products. 2020. <https://www.fda.gov/media/121751/download>
  - 18). European Network of GMO Laboratories. Definition of Minimum Performance Requirements for Analytical Methods of GMO Testing. 2015. [https://gmo-crl.jrc.ec.europa.eu/doc/MPR%20Report%20Application%2020\\_10\\_2015.pdf](https://gmo-crl.jrc.ec.europa.eu/doc/MPR%20Report%20Application%2020_10_2015.pdf)
  - 19). ISO. Biotechnology — Requirements for evaluating the performance of quantification methods for nucleic acid target sequences — qPCR and dPCR. ISO 20395:2019(en). <https://www.iso.org/obp/ui#iso:std:iso:20395:ed-1:v1:en>
  - 20). Sidstedt M, Rådström P, Hedman J. PCR inhibition in qPCR, dPCR and MPS-mechanisms and solutions. *Anal Bioanal Chem*. 2020;412(9):2009-2023.
  - 21). Schrader C, Schielke A, Ellerbroek L, Johne R. PCR inhibitors - occurrence, properties and removal. *J Appl Microbiol*. 2012;113(5):1014-26.
  - 22). Khan G, Kangro HO, Coates PJ, Heath RB. Inhibitory effects of urine on the polymerase chain reaction for cytomegalovirus DNA. *J Clin Pathol*. 1991;44(5):360-5.
  - 23). Demeke T, Jenkins GR. Influence of DNA extraction methods, PCR inhibitors and quantification methods on real-time PCR assay of biotechnology-derived traits. *Anal Bioanal Chem*. 2010;396(6):1977-90.
  - 24). Yedidag EN, Koffron AJ, Mueller KH, Kaplan B, Kaufman DB, Fryer JP, Stuart FP, Abecassis M. Acyclovir triphosphate inhibits the diagnostic polymerase chain reaction for cytomegalovirus. *Transplantation*. 1996;62(2):238-42.