

■ *gpt delta* マウス・ラットを用いる遺伝子突然変異試験に使用される大腸菌株

strain name	comment	reference
YG6020	For <i>gpt</i> assay (6-thioguanine selection), Kanamycin resistant	Nohmi et al. (1996) Environ.Mol.Mutagen., 28, 465-470.
<i>E. coli</i> C	For titer check (check of packaging efficiency)	Gossen et al. (1988) Nuc. Acids. Res. 16, 9343
WL95 (P2 lysogen)	For <i>Spi</i> <sup>-</sup> assay (confirmation of <i>Spi</i> <sup>-</sup> mutant phenotype)	Ikeda et al. (1995) Adv.Biophys., 31, 197-208.